

PRODUCT
CATALOG
2006

Finally a choice...

simple

reports of less pain

cost effective

versatile

Negative Pressure Wound Therapy



Versatile 1™ Wound Vacuum System*

The most cost effective and gentle wound vacuum on the market today!

The Versatile 1™ System Includes:

- The Versatile 1™ Pump with 15 ft Power Cord
- Small Disposable Canister (250 ml)
- Large Disposable Canister (800 cc)
- Pump Canister Connector
- Bacteria/Overflow Filters (5 ea)
- Canister Hoop Adapters
- Instructions for use

Item #100.0100

*Shown with optional Aqualert™



Versatile 1™ Accessories

(all sold individually)

VERSATILE 1™ PUMP STAND

This 5 leg, wide base stand for the Versatile 1™ enhances the convenience, safety and mobility of the Negative Pressure Wound Vacuum System. Tamper-proof screws secure the Versatile 1™ Wound Vacuum System to the stand. Also comes with a convenient basket that can hold canister and wound sealing kits. Two wheel locks keep the stand in place.

Item # 800.0162



CAR CONNECTOR ADAPTER

Make the Versatile 1™ mobile with a handy car connector for those long trips.

Item # 800.0101



VERSATILE 1™ CARRYALL

A handy carry all designed for the Versatile 1™, with additional pocket storage space.


Item # 800.0350



UNIVERSAL HOLDER

Slides into the back of the Versatile 1™. Can mount on bed rails, tables, IV poles, etc...

Item # 800.0205



Miller DermiVex™ Wound Drain

Introducing the Miller DermiVex™ Wound Drain. The principle of using low levels of negative pressure in wounds has been well established. Now a simple yet effective approach is available.

Less cutting.

Less fitting.

Item # 500.7026



pH Test Strips

15 strips per pack
25 packs per case

These pH Testing Strips provide an accurate way to measure the pH of a wound in vitro. Negative pressure wound therapy assessments may include pH testing as a quantitative way of measuring wound healing. pH range of 4.5-9.0

Item # 700.0806



Lanabalm™

Ultimate Care for Feet, Elbows, and Cracked Skin

LanaBalm™ is a rich, thick lotion of 100% medical grade lanolin that keeps moisture sealed in and restores the skin to its original condition. It is safe for use on sensitive skin and doctor recommended for diabetics.

4 oz tube

24 per case

Item # 700.0357



Aqualert™

In-line overflow alarm that can be used with all medical pumps in the BlueSky product line and fits most standard tubing. When fluid is detected in the tubing, the Aqualert™ activates an audible and visual alarm.

*Tubing & Bacteria/Overflow Protection Filter Included

Item # 700.0310



Small Disposable Canister

250 ml disposable canister with mechanical overflow valve.

Item # 700.9405

Small Hoop Adapter

Item # 800.5525


Wound Sealing Kit Components




Flat Drain – 10 mm




Round Drain – 10,19 FR



Channel Drain – 15 FR



Irrigation/Aspiration Drain



A. Non-Adherent Gauze

- Small 3" x 3"
- Medium 3" x 8"
- Large 3" x 8" (3 pack)

B. Anti-microbial Gauze

- Small 4" x 4"
- Medium 6" x 6" (2 pack)
- Large Roll

C. NPWT Transparent Adhesive Dressing

- Small 4" x 4.75"
- Medium 6" x 8"
- Large 8" x 12"

D. Complete Kit Components

- Sterile Towel/Drape
- Scissors
- No-Sting Skin Prep
- Sterile saline bullets
- Powder-free gloves

Wound Sealing Kit Components

Product Number	Wound Sealing Kit	Components *															
		Drains						Non-Adherent Gauze			Anti-Microbial Gauze			Transparent Adhesive Dressing			Tubing w/ Adaptor Connector
		10 mm Flat	10 FR Round	19 FR Round	15 FR Channel	Irrig/Asp.	DermiVex™	Small	Medium	Large	Small	Medium	Large	Small	Medium	Large	
500.2005	Chariker-Jeter® Basic	X															X
500.2045	Chariker-Jeter® Basic Channel				X												X
500.2060	Chariker-Jeter® Basic Round		X														X
500.3010	Chariker-Jeter® Basic Plus	X						X			X			X			X
500.3015	Chariker-Jeter® Basic Channel Plus				X			X			X			X			X
500.7777	Chariker-Jeter® Complete **	X						X			X		X		X	X	X
500.4010	Small Kit, 10 FR Round		X					X			X			2			
500.4011	Small Kit, Flat	X						X			X			2			
500.4050	Medium Kit, 15 FR Channel				X			X			X			2			
500.4051	Medium Kit, Flat	X						X			X			2			
500.4070	Large Kit, 19 FR Round			X					X			X			2		
500.4071	Large Kit, Flat	X							X			X			2		
500.4090	X-Large Kit, 19 Fr Round			X					X			2			3		
500.8743	Fill Kit							X			X			2			
500.1035	Wooding-Scott Basic					X											X
500.8888	Wooding-Scott Complete **					X			X		X		X		X	X	X
500.9999	Kremlin Kit	X															X
<i>Kremlin Kit also includes a Soft Suction Dome, Hard Suction Dome, stopper, and stretch garment</i>																	
500.7026	Miller DermiVex™ Drain					X											X

* All Kits include a drain clamp

** Complete kits also include scissors, gloves, sterile towel/drape, sterile saline, and No-Sting Skin Prep

Look for our Kit Selection Guide online soon!

Additional Kit Options

Miller DermiVex™

Simple to Apply

An easy to apply solution for wound therapy that provides negative pressure in the wound bed. By using the Miller DermiVex™ you can reduce the time, materials, and cost of wound dressings needed to fit a patient. The thin yet strong, clear silicone of the Miller DermiVex™ can be easily trimmed to size.


The convex drain is placed over the wound bed and negative pressure causes it to contract and seal the wound. No more difficult-to-apply adhesives! Once the wound is sealed, the Miller DermiVex™ works the same as any other wound drain by aspirating fluid from the wound while applying a comfortable and therapeutic level of negative pressure.

Features:

7 in x 4 in oval; shaped to fit most wounds

Clear, durable silicone for visual review

Functions as a drain and dressing



Miller DermiVex™ Wound Sealing Kit Contents:

1 ea—Silicone Wound Drain

1 ea—Tubing with Adapter Connector

Item # 500.7026



BACTERIA/OVERFLOW PROTECTION FILTER

The second level of overflow protection is the Bacteria/Overflow filter for the Versatile 1 Wound Vacuum System.

*Tubing included

Item # 700.6050



ODOR CONTROL FILTER

Assists with odor control. Must be used with the Bacteria/Overflow Filter. Not hydrophobic.

*Tubing included

Item # 700.0601



OSTOMY STRIP PASTE

10 Strips per Box

Item # 700.0850




Y CONNECTOR

Use the Y-Connector for multiple wound sites.

*Tubing included

Item # 700.2105



WOUND RULER

Disposable Ruler Packets

10 rulers per packet

100 packets per Case

Item # 700.0880



SOLIDIFIER

For use in canisters.

50 per Box

Item # 700.9081

Canister Kit Options



800 cc Disposable Canister




250 ml Disposable Canister



Canister with Pump/Canister Connector



Complete Canister Kit



Complete Small Canister Kit

Product Number	Canister Kits	# per Case	Components					
			800 ml Disposable Canister	250 ml Disposable Canister	Pump/Canister Connector	Antibacterial overflow filter	Tubing w/ Adaptor Connector	Solidifier
700.9005	Canister-Disposable	24	X					
700.9405	Disposable Small Canister	24		X				
700.9025	Canister, P/C connector, and anti-bacterial overflow filter	10	X		X	X		
700.9505	Canister, P/C connector, and anti-bacterial overflow filter PLUS Solidifier	10	X		X	X		X
700.9555	Complete Canister Kit*	10	X		X	X	X	
700.9565	Complete Canister Kit plus Solidifier*	10	X		X	X	X	X
700.9445	Completer Small Canister Kit*	10		X	X	X	X	

* Suggested for use with S,M,L,XL Wound Sealing Kits

The Versatile 1™ Wound Vacuum System


Negative Pressure Wound Therapy
Medicare Code E2402

Finally a cost-effective negative pressure wound therapy approach (www.npwt.com). Includes the Chariker-Jeter® Wound Sealing Kit designed by the pioneers in the field. Chariker-Jeter® kits last longer than conventional dressings and no tissue tearing or in-growth has been reported with these kits.

Simple
Cost Effective
Choices for Advanced Practitioners

BlueSky
Medical

BlueSky Medical
5924 Balfour Court, Suite 102
Carlsbad, CA 92008



PROVIDER INFORMATION HERE